




Diamond design



Back chamfer



transparent tubing
clear blood return

Disposable Injection Needles

Easy-Effective-Exactly

A Xinwell Member Company

AteTec

Intended use

Disposable endoscopy injection needle, passing through endoscopy channel into human body cavity, operates the injection of liquid after endoscopic mucosal biopsy and polypectomy.



Products Features

- The tube body is transparent, the blood return is clear, and the safety is guaranteed.
- The needle-point diamond design, double-edged, smooth needle entry.
- Good conductivity, the needle can be pulled out smoothly after the tube body is bent and circling.
- Secondary de-sheathing design, calmly deal with sticky needles.



Model	Diameter of outer tube(mm)	Length (mm)	Needle diameter (mm)	Needle length (mm)
ATE-ZSZ-18x1200x25x4	1.8	1200	25G	4
ATE-ZSZ-18x1800x25x6	1.8	1800	25G	6
ATE-ZSZ-23x1600x23x5	2.3	1600	23G	5
ATE-ZSZ-23x2300x21x4	2.3	2300	21G	4
ATE-ZSZ-23x1800x21x5	2.3	1800	21G	5
ATE-ZSZ-23x1800x21x6	2.3	1800	21G	6
ATE-ZSZ-23x1800x23x4	2.3	1800	23G	4
ATE-ZSZ-23x1800x23x5	2.3	1800	23G	5
ATE-ZSZ-23x1800x23x6	2.3	1800	23G	6
ATE-ZSZ-23x1800x25x4	2.3	1800	25G	4
ATE-ZSZ-23x1800x25x5	2.3	1800	25G	5
ATE-ZSZ-23x1800x25x6	2.3	1800	25G	6
ATE-ZSZ-23x2300x21x4	2.3	2300	21G	4
ATE-ZSZ-23x2300x21x5	2.3	2300	21G	5
ATE-ZSZ-23x2300x21x6	2.3	2300	21G	6
ATE-ZSZ-23x2300x23x4	2.3	2300	23G	4
ATE-ZSZ-23x2300x23x5	2.3	2300	23G	5
ATE-ZSZ-23x2300x23x6	2.3	2300	23G	6
ATE-ZSZ-23x2300x25x4	2.3	2300	25G	4
ATE-ZSZ-23x2300x25x5	2.3	2300	25G	5
ATE-ZSZ-23x2300x25x6	2.3	2300	25G	6

JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone
213000 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

+86-519 8528 8988

+86-519 8528 8060

A Xinwell Member Company

AteTec