



Rotatable/non-rotatable

Oval

Hexagon

Crescent



Disposable Polypectomy Snares
Easy-Effective-Exactly

A Xinwell Member Company

AteTec

Intended use

The device is used endoscopically in the removal of polyps, and tissue from within the GI tract.

Products Features

- The control handle can make the snare rotate synchronously and effectively improve the success rate of the operation.
- Using imported special high-strength medical stainless steel wire material, cutting faster, Smoother and more uniform, effectively reducing bleeding damage.
- The specifications and models are rich and varied to meet the various clinical needs of doctors.
- It is compatible with all mainstream high-frequency electrical devices on the market.

Model	Shape of Loop	Diameter of outer tube(mm)	Length (mm)	Loop width (mm)
ATE-QTQ-TY-23 × (1600/1800/2300) × 10	Oval	2.3	(1600/1800/2300)	10
ATE-QTQ-TY-23 × (1600/1800/2300) × 15	Oval	2.3	(1600/1800/2300)	15
ATE-QTQ-TY-23 × (1600/1800/2300) × 20	Oval	2.3	(1600/1800/2300)	20
ATE-QTQ-TY-23 × (1600/1800/2300) × 25	Oval	2.3	(1600/1800/2300)	25
ATE-QTQ-TY-23 × (1600/1800/2300) × 30	Oval	2.3	(1600/1800/2300)	30
ATE-QTQ-TY-23 × (1600/1800/2300) × 35	Oval	2.3	(1600/1800/2300)	35
ATE-QTQ-LJ-23 × (1600/1800/2300) × 10	Hexagon	2.3	(1600/1800/2300)	10
ATE-QTQ-LJ-23 × (1600/1800/2300) × 15	Hexagon	2.3	(1600/1800/2300)	15
ATE-QTQ-LJ-23 × (1600/1800/2300) × 20	Hexagon	2.3	(1600/1800/2300)	20
ATE-QTQ-LJ-23 × (1600/1800/2300) × 25	Hexagon	2.3	(1600/1800/2300)	25
ATE-QTQ-LJ-23 × (1600/1800/2300) × 30	Hexagon	2.3	(1600/1800/2300)	30
ATE-QTQ-LJ-23 × (1600/1800/2300) × 35	Hexagon	2.3	(1600/1800/2300)	35
ATE-QTQ-BY-23 × (1600/1800/2300) × 10	Crescent	2.3	(1600/1800/2300)	10
ATE-QTQ-BY-23 × (1600/1800/2300) × 15	Crescent	2.3	(1600/1800/2300)	15
ATE-QTQ-BY-23 × (1600/1800/2300) × 20	Crescent	2.3	(1600/1800/2300)	20
ATE-QTQ-BY-23 × (1600/1800/2300) × 25	Crescent	2.3	(1600/1800/2300)	25
ATE-QTQ-BY-23 × (1600/1800/2300) × 30	Crescent	2.3	(1600/1800/2300)	30
ATE-QTQ-BY-23 × (1600/1800/2300) × 35	Crescent	2.3	(1600/1800/2300)	35
ATE-QTQ-TY-18 × (1600/1800) × 10/15	Oval	1.8	(1600/1800)	10/15
ATE-QTQ-LJ-18 × (1600/1800) × 10/15	Hexagon	1.8	(1600/1800)	10/15
ATE-QTQ-BY-18 × (1600/1800) × 10/15	Crescent	1.8	(1600/1800)	10/15

JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone
213000 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

 +86-519 8528 8988

 +86-519 8528 8060

A Xinwell Member Company

AteTec