



Disposable Hemostatic Clips

Easy-Effective-Exactly

A Xinwell Member Company

AteTec

Intended use

This hemostatic clips are compatible with endoscopes ,clamps are placed in the upper and lower gastrointestinal tracts to stop bleeding caused by various reasons and wounds after endoscopic surgery.

Products Features

- Ergonomically designed handle is more accurate and convenient to operate.
- Smooth coating design to avoid damage to the instrument channel of the endoscope .
- The handle can rotate synchronously with the clamp device, which can meet the needs of different lesions on the release angle and direction of the clamp device in surgery.
- Variety spans of 10mm、 13mm and 16mm,which fully meet the different requirements. Excellent clamping strength,safety and reliable.
- The clip of the clip device adopts the material with excellent toughness, which can reduce the injury of the affected part during the operation. There will be a prompt sound after release, which can improve the safety performance of the instrument and reduce the risk in the operation.
- Two types are available (with or without infusion port).

Model	Diameter of Outer Tube (mm)	Clip Opening (mm)	Working length (mm)	Characteristic
ATE-ZXJ-B-2.3×1950×135×9	2.3	9	1950	With Infusion Port
ATE-ZXJ-B-2.3×2350×135×9	2.3	9	2350	With Infusion Port
ATE-ZXJ-B-2.3×1950×135×11	2.3	11	1950	With Infusion Port
ATE-ZXJ-B-2.3×2350×135×11	2.3	11	2350	With Infusion Port
ATE-ZXJ-B-2.3×1950×135×13	2.3	13	1950	With Infusion Port
ATE-ZXJ-B-2.3×2350×135×13	2.3	13	2350	With Infusion Port
ATE-ZXJ-B-2.3×1950×135×16	2.3	16	1950	With Infusion Port
ATE-ZXJ-B-2.3×2350×135×16	2.3	16	2350	With Infusion Port
ATE-ZXJ-A-2.3×1950×135×9	2.3	9	1950	N/A
ATE-ZXJ-A-2.3×2350×135×9	2.3	9	2350	N/A
ATE-ZXJ-A-2.3×1950×135×11	2.3	11	1950	N/A
ATE-ZXJ-A-2.3×2350×135×11	2.3	11	2350	N/A
ATE-ZXJ-A-2.3×1950×135×13	2.3	13	1950	N/A
ATE-ZXJ-A-2.3×2350×135×13	2.3	13	2350	N/A
ATE-ZXJ-A-2.3×1950×135×16	2.3	16	1950	N/A
ATE-ZXJ-A-2.3×2350×135×16	2.3	16	2350	N/A

JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone
213000 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

 +86-519 8528 8988

 +86-519 8528 8060

A Xinwell Member Company

AteTec