



Disposable Bite Blocks

A Xinwell Member Company

AteTec

Intended use

Disposable Bite Blocks designed to help ensure that the mouth of the patient remains open during endoscopic procedures. Provides protection for patients and the endoscopic equipment.

Products Features

- Different models and colors are available.
- Smooth surface to enhance patient comfort.
- Standard mouthguards with strap suitable for adult to improve safety and convenience.

Picture	Model	Width (mm)
	ATE-YO-02	45±5
	ATE-YO-03	26±5

JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone
213000 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

 +86-519 8528 8988

 +86-519 8528 8060

A Xinwell Member Company

AteTec