



Disposable Endoscope Valves Sets

A Xinwell Member Company

AteTec

Intended use

Disposable Endoscope Valves Sets are used for multiple endoscopes working channels or ports to control the flow of liquids, gases, and other substances.

- **Disposable Air/Water Valves** are intended to control the inflow of gases, water and prevent backflow.
- **Disposable Suction Valves** are intended to enable suction and prevent inflow of air.
- **Disposable Biopsy Valves** are intended to be used at the valve opening of the endoscopic biopsy forceps to prevent gas and fluid leakage during endoscopic surgery.

Products Features

- **Single-use**
Aviod cross contamination
- **Ethylene oxide sterilization**
Safety is guaranteed
- **Easy assemble**
Not easy to fall off
- **Compatible with various series of endoscopes**



Specifications

	Included product	Included product model	QTY	Compatible endoscopes models
ATE-FTZ-A	Disposable Air/Water Valve	ATE-SQF-00	1	Olympus 140/160/180/190/240/260/290 Series of endoscopes
	Disposable Suction Valves	ATE-XYF-00	1	
	Disposable Biopsy Valves	ATE-HJF-00	1	
ATE-FTZ-B	Disposable Air/Water Valve	ATE-SQF-00	1	
	Disposable Suction Valves	ATE-XYF-00	1	
ATE-HJF-00	Disposable Biopsy Valves	/	1	
ATE-HJF-01	Disposable Biopsy Valves	/	1	Pentax series endoscopes

▲ If there is any change in product specification, design and accessories, please refer to the product registration information.