



Disposable Grasping Forceps

Visualization Organization Efficiency

A Xinwell Member Company

AteTec

Intended use

The disposable grasping forceps is used in conjunction with the soft endoscopy, entering into the human body cavity such as respiratory tract, esophagus, stomach, intestine and so forth through the endoscopy channel, to take tissues, stones, foreign bodies and stent out.

Products Features

- Integral rotation handle.
- With good flexibility, accurately and firmly grasp the foreign body.
- ESD and EMR's indispensable partner, the clots, polyps, tissues can be thoroughly retrieved.
- Flexible clamp head in opening and closing.
- Strength of product connection over 50N.



	Model	Diameter of the outer tube(mm)	Length(mm)
	Rat teeth with alligator		
	ATE-YWQ-CXDZ-18×1600	1.8	1600
	ATE-YWQ-CXDZ-23×1600	2.3	1600
	ATE-YWQ-CXDZ-23×2300	2.3	2300
	3 Prongs		
	ATE-YWQ-3ZX-18×1600	1.8	1600
	ATE-YWQ-3ZX-23×1600	2.3	1600
	ATE-YWQ-3ZX-23×2300	2.3	2300
	4 Prongs		
	ATE-YWQ-4ZX-18×1600	1.8	1600
	ATE-YWQ-4ZX-23×1600	2.3	1600
	ATE-YWQ-4ZX-23×2300	2.3	2300
	5 Prongs		
	ATE-YWQ-5ZX-18×1600	1.8	1600
	ATE-YWQ-5ZX-23×1600	2.3	1600
	ATE-YWQ-5ZX-23×2300	2.3	2300

JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone
213000 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

+86-519 8528 8988

+86-519 8528 8060

A Xinwell Member Company

AteTec