



Disposable Multi-band Ligator


Esophageal and gastric varices ligation and
internal hemorrhoid ligation

A Xinwell Member Company

AteTec


Intended use

It is used in endoscopes for esophageal and gastric varices ligation and internal hemorrhoid ligation.



Products Features

- Wide surgical field and high definition are convenient for doctors to quickly identify lesions.
- Ergonomic handle design, comfortable and convenient operation, strong controllability.
- Good elasticity and resilience of ligature ring improves the reliability of the ligature.
- Firm structure and easy operation, ensuring the smooth release of the ligature ring.



Model	Working Length (mm)	Number of Ligature Ring	Outer Diameter of Endoscope Tip (mm)
ATE-TZQ-1530×4×B	1530	4	Φ10 — Φ11.9
ATE-TZQ-1530×5×B	1530	5	Φ10 — Φ11.9
ATE-TZQ-1530×6×B	1530	6	Φ10 — Φ11.9
ATE-TZQ-1530×7×B	1530	7	Φ10 — Φ11.9
ATE-TZQ-1530×8×B	1530	8	Φ10 — Φ11.9
ATE-TZQ-1530×9×B	1530	9	Φ10 — Φ11.9

JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone
213000 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

 +86-519 8528 8988

 +86-519 8528 8060

A Xinwell Member Company

AteTec